
Videolaryngoskop

VX-30 & VX-35

EN	Video laryngoscope VX-30 & VX-35
FR	Vidéolaryngoscope VX-30 & VX-35
IT	Videolaringoscopio VX-30 e VX-35
ES	Videlaringoscopio VX-30 y VX-35
CS	Videolaryngoskop VX-30 a VX-35
NL	Videolaryngoscoop VX-30 & VX-35
SV	Videolaryngoskop VX-30 & VX-35
PL	Laryngoskop Wideo VX-30 & VX-35
PT	Videolaringoscópio VX-30 e VX-35

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1. USB-C (Laden & Externe Wiedergabe)
2. Display
3. Power-Button
4. Foto-Video-Button
5. Handgriff
6. Arretierung Einwegspatel
7. HD-Kamera

EN – Product Overview (S. 22)	FR - Présentation du produit (S. 39)	IT – Panoramica del prodotto (S.57)
1. USB-C (Charging & Streaming) 2. Display 3. Power-Button 4. Photo-Video-Button 5. Handle 6. Disposable Blade Lock 7. HD-Camera	1. USB-C (chargement et lecture externe) 2. Écran 3. Bouton d'alimentation 4. Bouton photo/vidéo 5. Poignée 6. Verrouillage de la spatule jetable 7. Caméra HD	1. USB-C (ricarica e streaming) 2. Display 3. Pulsante di accensione 4. Pulsante foto-video 5. Impugnatura 6. Blocco lama monouso 7. Fotocamera HD

ES – Descripción general (S. 75)	CS – Přehled produktu (S. 93)	NL – Productoverzicht (S. 111)
1. USB-C (carga y transmisión) 2. Pantalla 3. Botón de encendido 4. Botón de foto y video 5. Asa 6. Bloqueo de cuchilla desechable 7. Cámara HD	1. USB-C (nabíjení a streamování) 2. Displej 3. Tlačítko napájení 4. Tlačítko foto/video 5. Rukojeť 6. Zámek jednorázové čepele 7. HD kamera	1. USB-C (opladen en streamen) 2. Display 3. Aan / uit-knop 4. Foto-/videoknop 5. Handgreep 6. Vergrendeling voor wegwerpmesjes 7. HD-camera

SV – Produktöversikt (S. 129)	PL – Przegląd produktu (S. 147)	PT – Visão geral do produto (S. 165)
1. USB-C (laddning och streaming) 2. Display 3. Strömröret 4. Foto-/videoknapp 5. Handtag 6. Lås för engångsblad 7. HD-kamera	1. USB-C 2. Wyświetlacz 3. Przycisk zasilania 4. Przycisk zdjęć i filmów 5. Uchwyt 6. Blokada jednorazowego ostrza 7. Kamera HD	1. USB-C (laddning och streaming) 2. Ecrã 3. Botão liga/desliga 4. Botão foto/vídeo 5. Pega 6. Trava da lâmina descartável 7. Câmara HD

DE - VIDEO LARYNGOSKOP VX-30 & VX-35

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EN - VIDEO LARYNGOSCOPE VX-30 & VX-35

This product is intended for use by licensed physicians and trained medical personnel only.

Please read this instruction manual carefully before using this product and keep it in a safe place for future reference.

During use, the operating instructions must be strictly followed, and maintenance and care must be conducted properly.

Paragraphs marked "Warning" must be read carefully and observed to prevent damage to the device and injury to the operator or patient.

If problems arise during operation of the device, please contact your local distributor or our company directly. We are committed to providing you with first-class service and support.

Commitment: Should troubleshooting be necessary, we are committed to providing the user with all necessary technical documentation containing detailed information upon request.

Note: We improve our Product continuously. The actual product may differ from this manual in terms of specifications or appearance.

Scope of application

The VX-30 & VX-35 video laryngoscopes are used in clinical and preclinical applications to lift the patient's epiglottis, exposing the glottis and enabling medical professionals to perform precise intubation. They can also be used for intraoral examinations and treatments.

Function

The video laryngoscope is a medical device designed to assist with intubation. Its rigid blade stabilizes the tongue and allows the user to visualize the laryngeal structures.

1. 3.2 Model and Specifications

Table 1; Display


Type	Model	Display size	Touchscreen	REF
Display	Monitor MX-30	3 inches	X	I-VL-02-0010
	Monitor MX-35	3.5 inches	✓	I-VL-02-0011

Type	Model	Type	REF
Handle	HX-M handle	MAC	I-VL-02-0012

The optional accessories for the video laryngoscope, along with their specifications and model numbers, are listed in Table 2.

Table 2; Disposable blades

Name	Type	Dimensions				REF
		L $\pm$ 10%	B1 $\pm$ 10%	B2 $\pm$ 10%	H $\pm$ 10%	
Disposable Blade	MAC 4	132	18	29	52	I-VL-03-0013
	DA 3	116	16	24	51	I-VL-03-0012
	MAC 3	119	16	26	47	I-VL-03-0011
	MAC 2	109	13	22	43	I-VL-03-0010
	MAC 1	104	12	20	39	I-VL-03-0009
	INF 1	106	10	16	41	I-VL-03-0008
	INF 0	102	9	15	38	I-VL-03-0007

- This device is only compatible with VisionHEART disposable blades.
- Material: Polycarbonate (PC)
- Packaging: Sterile packaging, for single use only. Device and packaging are latex-free.

1.4 Compliance

The device complies with the safety standards set out in ISO 7376, EN 60601-1, and EN 60601-1-2. It bears the CE mark and meets the requirements of European Commission Directives 93/42/EEC and 2007/47/EC on medical devices. In addition to the above standards, the sterility of VisionHEART blades has been confirmed by the TÜV certification body.

In accordance with section 8.3 of EN 60601-1 on touchable parts, VisionHEART disposable blades are classified as type BF.

Users should contact their VisionHEART representative or distributor to ensure safe and responsible disposal in accordance with applicable regulations on waste electrical and electronic equipment (WEEE) and, where applicable, future environmental regulations (including, but not limited to, Directive 2011/65/EU, implemented in the respective country or region of use).

The device complies with the safety standards set forth in ISO 7376, EN 60601-1, and EN 60601-1-2. It bears the CE mark and meets the requirements of European Commission Directives 93/42/EEC and 2007/47/EC on medical devices. In addition to the above standards, the sterility of VisionHEART blades have been confirmed by the TÜV certification body.

In accordance with Section 8.3 of EN 60601-1 on touchable parts, VisionHEART disposable blades are classified as Type BF.

Users should contact their VisionHEART representative or distributor to ensure safe and responsible disposal in accordance with applicable regulations on waste electrical and electronic equipment (WEEE) and, where applicable, future environmental regulations (including, but not limited to, Directive 2011/65/EU, implemented in the respective country or region of use).

1.5 Environmental Conditions

Application conditions	Ambient temperature	-10~+50°C
	Relative humidity	≤80
	Atmospheric pressure	700 hPa to 1060 hPa
Storage and transport conditions	Ambient temperature	-20°C~+55°C
	Relative humidity	≤93%
	Atmospheric pressure	500 hPa to 1060 hPa

2.0 Functionality

2.1 Check before use

Open the packaging and check for transport damage. If there are any signs of damage, do not use the device and inform your VisionHEART representative or distributor.

The device is non-sterile. The disposable blades are in a sterile condition. Ensure that handling procedures comply with local, national, and harmonized standards.

2.2 Battery

Before use, ensure that the battery is sufficiently charged. If not used for an extended period, recharge every 2-3 months.

2.3 Mounting the monitor and locking the disposable blade

Insert the monitor into the handle until it clicks into place. The lock holds the monitor securely in the handle. Before inserting the disposable blade, open the sterile packaging and insert the blade into the handle. After the locking mechanism has **audibly** clicked into place, carefully move the blade in all directions to ensure that it is secure and will not loosen or come out.

VisionHEART disposable blades are intended for single use only and must be discarded after each use on a patient.


Ensure that local handling regulations are always followed when handling the sterilized device.


VisionHEART disposable blades are only compatible with the VisionHEART brand video laryngoscope.

Do not use the disposable blade if the packaging shows signs of damage.

Press the monitor into the handle until it clicks into place.

Slide the disposable blade onto the handle until...


... it audibly clicks into place.


2.4 Using the video laryngoscope

Place the patient in the optimal "sniffing position"

Carefully insert the blade into the oral cavity


Insert the device along the midline while moving the tongue to the left.


Advance the tip of the blade to the epiglottis.


When the epiglottis is visible, position it so that the glottis is visible in the center

The screen allows real-time visualization of the intubation process.


The image on the screen can be used to check the correct insertion depth of the endotracheal tube.

Insert the tube slowly and monitor its progress either orally or on the screen.

2.5 Clinical use


Insert the mandrin/stylet into the tube and shape it according to the curvature of the disposable blade.


Position the patient's head in the optimal "sniffing position" for intubation


Rotate the blade along the contour of the tongue until the epiglottis becomes visible in the center of the screen.


Carefully insert the blade into the oral cavity.


Carefully rotate the laryngoscope until the epiglottis is lifted.


When optimally positioned, the glottis is visible in the upper center of the screen


Insert the tube from the right side of the patient's mouth and gently advance it along the intubation guide of the blade.

If the intubation is performed correctly, you will see the tube moving forward along the right guide until its tip reaches the front of the glottis.

Hold the mandrin securely and rotate the tube into the glottis. The mandrin should not enter the glottis. Once the tube has passed the vocal cords, remove the mandrin.

The image on the screen provides an accurate visualization of the insertion depth of the tube.

3. Maintenance and sterilization

The device should be decontaminated immediately after each use on a patient. When handling and cleaning contaminated items, hospital guidelines and protocols must be followed.

The following cleaning methods comply with globally recognized standards for medical facilities, developed in accordance with relevant legal requirements and taking into account AAMI TIR12 (Design, Testing, and Labeling of Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers).

The instructions below have been certified and confirm their suitability for preparing medical devices for reuse. It is the responsibility of the reprocessing or to ensure that the reprocessing procedures, including the actual use of equipment and materials and the personnel in the reprocessing facility, achieve the intended results. This typically requires validation of the relevant procedures and routine monitoring.

The following two methods are approved by the manufacturer:

Product	Cleaning	Disinfection
	Delivered sterile & single use only!	Delivered sterile & single use only!
	Manual cleaning – without immersion Protocol for cleaning with enzymatic cleaning agent	Manual cleaning – without immersion Disinfectant wipe with 70% ethanol swab

Avoid using strong acids, highly oxidizing disinfectants, or disinfection methods that exceed the specified environmental tolerances of the product.

Ensure that the device is completely dry before use.

3.1 Cleaning

Follow the manufacturer's instructions for preparing or diluting the enzymatic cleaning solution. Disassemble the handle and blade (including the camera), use purified water to prepare the solution, and thoroughly wipe all exposed surfaces of the video laryngoscope. Both the handle and blade (including the camera) must be wiped for at least 2 minutes each.

Allow the video laryngoscope to air dry, then wipe all exposed surfaces of the video laryngoscope with sterile cotton balls until no residue remains.

3.2 Disinfection

1. Disinfect the video laryngoscope with disinfectant solution and allow it to air dry for 20 minutes.
2. Wipe all exposed surfaces of the video laryngoscope with sterile gauze soaked in sterile water.
3. Clean the handle and handpiece (including camera) for 2 minutes each with 70% ethanol swabs.

4. Storage and transport

The following precautions must be observed when storing and transporting the device:

1. Ensure that contaminated devices are stored separately in a suitable container in accordance with local hospital regulations. The container used to transport the device should be durable and suitable for decontamination.
2. Maintain cleanliness during transport and avoid contamination; be sure to protect the device from rain and snow.
3. Store the device in a well-ventilated indoor area that is free of harmful gases and external biological contaminants.

5. Troubleshooting

Malfunction	Cause	Troubleshooting
Poor image quality	The camera is dirty.	Clean the camera
The device cannot be turned on	Battery empty	Charge the battery
Device is turned on, but no image appears on the screen	Poor contact between handle and display	Clean the contact points between the handle and the display and restart the device.

If the above steps do not improve the image quality, contact your VisionHEART representative or dealer for diagnosis and repair.

6. Warranty

The manufacturer of the VisionHEART device provides a 24-month warranty.

Important note: To prevent the transmission of infection and ensure the safety of all employees who manage medical devices, all devices must be thoroughly cleaned and subjected to high-level disinfection before being returned to our company for repair. If the device has been used on HBsAg-positive patients or other cases of infection, you must expressly inform the manufacturer/distributor.

Note:

Before each use, check the video laryngoscope and disposable blades for rough surfaces, sharp edges, or protrusions that could pose a safety risk when inserted into the human body. This product must not be used simultaneously with other electrical medical devices.

Manufacturer's declaration

The "VisionHEART" video laryngoscope (this device) complies with the IEC60601-1-2:2012 standard ["Medical electrical equipment – Part 1-2: General requirements for safety including essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests"].

In addition, the sections of IEC60601-2-18:2009 ["Medical electrical equipment – Part 2-18: Particular requirements for the safety including essential performance of endoscopic equipment"] also apply to this device insofar as they refer to IEC60601-1-2:2020.

The following table contains the manufacturer's declaration regarding electromagnetic emissions and electromagnetic immunity.

Table 1 – Electromagnetic emissions
This device is intended for use in the electromagnetic environment specified below. The user should ensure that the device is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
HF Emission GB4824	Group 1	The video laryngoscope is suitable for use in all facilities that are not domestic and are not directly connected to the public low-voltage power supply for residential buildings. The video laryngoscope uses high-frequency energy exclusively for its internal functions. Therefore, its RF emissions are very low, and the likelihood of it causing interference with nearby electronic devices is minimal.
HF Emission GB4824	Class A	
Harmonic emission GB 17625.1	Not applicable	
Voltage fluctuations / Flicker emission GB 17625.2	Not applicable	

Table 2 – Electromagnetic immunity			
This device is intended for use in the electromagnetic environment specified below. The user should ensure that the device is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance	Electromagnetic environment - Guidance
Electrostatic discharge GB/T 17626.2	±6kV contact discharge, ±8kV air discharge	±6kV contact discharge, ±8kV air discharge	The floor should be made of wood, concrete, or ceramic tiles. If the floor is covered with synthetic materials, the relative humidity should be at least 30%.

Fast electrical transients/bursts GB/T17626.4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	This clause is not applicable because the device is battery-powered.
Surge voltage GB/T 17626.5	± 1 kV conductor to conductor (L-L) ± 2 kV conductor to ground (L-G)	Not applicable	This clause is not applicable because the device is battery-powered.
Voltage dips, short-term interruptions, and voltage changes on the input power supply lines GB/T 17626.11	$< 5\% U_T$, duration 0.5 cycles ($> 95\%$ dip at U_T) $40\% U_T$, duration 5 cycles (60% dip at U_T) $70\% U_T$, duration 25 cycles (30% drop in U_T) $< 5\% U_T$, duration 5 seconds ($> 95\%$ drop in U_T)	Not applicable	This clause is not applicable because the device is battery-powered.
Magnetic field of the mains frequency (50Hz/60Hz) GB/T 17626.8	3A/m	3A/m	The power frequency magnetic field should have the typical characteristics and levels found in commercial or clinical environments.
Note: U_T refers to the AC voltage of the power grid before applying the test voltage.			

Manufacturing license no.

JS FDA Medical Device Manufacturing License 20190027

Medical device registration number

Su Xie Zhu Zhun 20212081332

Manufacturing license no.

JS FDA Medical Device Production License 20190027

Technical requirement number for products

Su Xie Zhu Zhun 20212081332

Name of registrant/manufacture	Jiangsu Yongle Medical Technology Co., Ltd.
Address of registrant/manufacture	No. 1 Jiaodang Industrial Cluster, Zhangqiao Town, Taixing City, Jiangsu Province, China
Contact details of the registrant/manufacture	+86-523-86550008
Postal code	225400
Name of customer service center:	Jiangsu Yongle Medical Technology Co., Ltd
Address	No. 1 Jiaodang Industrial Cluster, Zhangqiao Town, Taixing City, Jiangsu Province, China
Phone	0523-86550008 / Fax: 0523-86550006
Date of revision of the manual	August 1, 2024

Information card on hazardous substances in electronic and electrical products

1. Label for controlling environmental pollution from electronic and electrical products:

The following label indicates that this device does not contain any substances whose content exceeds the permissible limits.

2. Table of hazardous substances in electronic and electrical products: (Only for orange EFUP labeling)
No table of hazardous substances is required for this device. The device does not contain any substances whose content exceeds the permissible limits.

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VisionHEART GmbH
Everbiz GmbH
Jiangsu Yongle Medical Technology Co. Ltd.