



VIDEO LARYNGO SCOPE

BF-DP

versatile applications.

Various disposable blades for MAC, MIL, and Difficult Airway.

Simple documentation.

Capture images and videos with one-button operation.

User-friendly operation.

Anti-Fog Camera, mini HDMI output for training.



long-lasting Battery

Integrated battery with more than 120-minute runtime. Rechargeable.

training conditions

Mini-HDMI for live streaming to an external monitor and USB-C for file transfer to a computer or laptop. Perfect for training purposes

simple documentation

Document the intubation with the integrated memory card, which stores over 400 video recordings or more than 3 million images.

versatile applications

Great flexibility with a wide range of disposable blades for different patient groups. MIL 0, MAC 1, MAC 2, MAC 3, and MAC 4, as well as a specialized blade for difficult intubations (Difficult Airway).

user friendly operation

Easy to use with its integrated touchscreen. Adjust illumination, contrast and color management

designed for right- and left-handed users

Enabled by clear and central arrangement of the screen

high-quality materials

Monitor and handle are made from a special aluminum alloy, ensuring a high-quality feel and easy cleaning and disinfection.

HD camera

With excellent resolution of 1280 x 720 px



available configurations.

the video laryngoscope BF-DP is available in two different sets: one for adults (MAC) and one for neonates (MIL). Handles are also available separately, allowing users to complete the other set and cover both patient groups.



MAC Version



MIL Version

components & consumables.



3.5" Touchscreen Monitor



MAC Handle for MAC 1 -
MAC 4 & DMA



Neonate Handle for
MILO blade



MIL 0

MAC 1

MAC 2

MAC 3

MAC 4

Difficult
Airway

technical specification.

Parts	Description	Technical Index
Monitor	Screen	3.5" Touchscreen
	Resolution	960 x 640 (RGB)
	Viewing angle	> 70°
	Power	< 3W
	Rotation angle Left/right Front/Back	270° - 270° 190 - 190°
	Connectivity	Mini HDMI simultaneously
Camera	Illuminance	> 150 LUX
	Resolution	1280 x 720
Blade	Waterproof	IPX7
Image / Video Function	Image / Video function	Yes (with playback)
	Output	USB-C to Computer
	Memory Type	MicroSD 128GB
Battery	Battery Type	Recharchable Lithium Battery
	Capacity	3.7V / 1450 mAh
	Lifecycle of battery	approx. 500 times
	Battery working time	> 120 minutes
	Charging time	< 2 hours
	Port	USB-C

REF	Description	Application	
I-VL-01-0001	SET MAC	Monitor + Handle MAC	SET
I-VL-01-0002	SET MIL	Monitor + Handle MIL	
I-VL-02-0001	3.5" Touch- screen	LCD Monitor	COMPONENT
I-VL-02-0002	Adult Handle for MAC 1- 4 blades & DMA		
I-VL-02-0003	Neonate Handle for MIL 0 blade		
I-VL-03-001	Miller 0	Neonate	CONSUMABLES
I-VL-03-002	MAC 1	Pediatric	
I-VL-03-003	MAC 2	Small Adult	
I-VL-03-004	MAC 3	Adult	
I-VL-03-005	MAC 4	Large Adult	
I-VL-03-006	DMA	Difficult Airway	

Get in touch with you local Distributor: