

VIDEOLARYNGOSCOPE

VX-30 & VX-35



Product overview.

REF			Description	
I-VL-01-0003	Starter-Kit Videolaryngoscope	VX-30	Monitor 3 Inch + Handle + Bag	SET
I-VL-01-0004	Starter-Kit Videolaryngoscope	VX-35	Monitor 3.5 Inch + Handle + Bag	
I-VL-02-0010	Videolaryngoscope Monitor	MX-30	Monitor 3 Inch	COMPONENTS
I-VL-02-0011	Videolaryngoscope Monitor	MX-35	Monitor 3.5 Inch	
I-VL-02-0012	Videolaryngoscope Handle	HX-M	Handle	
I-VL-04-0006	Bag for Videolaryngoscope	VX-30 & VX-35		
I-VL-03-0007	Single-Use Laryngoscope Blade	INF 0	Neonates	CONSUMABLES
I-VL-03-0008	Single-Use Laryngoscope Blade	INF 1	Infants	
I-VL-03-0009	Single-Use Laryngoscope Blade	MAC 1	Pediatric	
I-VL-03-0010	Single-Use Laryngoscope Blade	MAC 2	Adult, small	
I-VL-03-0011	Single-Use Laryngoscope Blade	MAC 3	Adult	
I-VL-03-0012	Single-Use Laryngoscope Blade	MAC 4	Adult, large	
I-VL-03-0013	Single-Use Laryngoscope Blade	MAC DA3	Difficult Airway	

Rechargeable battery

Integrated and rechargeable battery with up to 180 minutes of operating time.

Plug & Play touchscreen

Intuitive operation via the touchscreen monitor (only 3.5") - available with a 3" or 3.5" display for maximum user comfort.

Documentation & access control

Document the intubation on the integrated memory card (over 400 videos / 3 million images). With the optional password lock, you can reliably secure your case records against unauthorised access.

Protected connector

Designed for user comfort - the recessed connection between blade and handle eliminates pressure points during intubation.

Made for all patients

Highly flexible due to disposable blades for different patient groups: INF 0, INF 1, MAC 1, MAC 2, MAC 3, MAC 4 and MAC DA3 – especially for difficult intubations (difficult airway).

Preformed guide channel

Intuitive and controlled tube guidance due to a guided channel integrated into the disposable blades.

Optimized for training purposes

Available for MX-35: HDMI (USB-C) for live streaming to an external monitor and USB-C for file transfer to a computer or laptop.

Ergonomic handle

Optimized for both left- and right-handed users – the monitor always remains centrally aligned, ensuring maximum visibility and ease of operation.

HD camera

2 MP CMOS camera with a resolution of 1280 x 720 px.



Reliable. Protected. Hygienic.

With its IP67 protection rating, the VX Series video laryngoscope is reliably protected in everyday clinical use as well as under extreme operating conditions. It is fully resistant to dust, liquids, and disinfectants.

The robust, fully sealed design minimizes the risk of cross-contamination and supports the highest hygiene standards. This ensures that the device not only remains fully functional, but also actively contributes to patient safety – wherever medicine leaves no room for failure.

Plug & Play Monitor.



The VX Series video laryngoscopes feature an intuitive plug-and-play interface.

Attach. Power on. Intubate.

Both monitors share the same intuitive user interface, outstanding image quality, and password-protected documentation functions.

The 3.5" monitor (MX-35) additionally includes an HDMI (USB-C) connection – providing even greater flexibility for training and everyday clinical use.

Consumables.

Preformed guide channel

Maximum control and safety: The preformed guiding channel directs the tube precisely through the camera's field of view, reduces unnecessary movements, and ensures a smooth, atraumatic intubation process.



The „curved Miller Blade“

Optimized tube guidance and seamless compatibility – the INF 0 and INF 1 blades are based on the classic MIL 0 & 1 design but have been refined with a curved geometry that perfectly complements the ergonomically curved handle.



STERILE

3 Years



Single-use